**Polyflush Syringe SS**  
0,9% sodium chloride (0,9% NaCl)  
Single sterile (SS) Sterile Filter (S) Path  
INSTRUCTIONS FOR USE

**MATERIALS USED**  
Polypropylene, SBR, Sodium Chloride

**INDICATION AND CLINICAL USE**  
Polyflush syringe with Saline solution is intended to be used only for flushing of indwelling vascular access devices (VAD's) for maintaining patency of VAD's.

**CONTRAINDICATIONS**  
None Known

**ADVERSE REACTION**  
Adverse reactions are known when the product is indicated.

**DIRECTIONS FOR USE**  
Use aseptic technique throughout the procedure.

1. Visually inspect product for particulate matter and discoloration prior to use according to Warning and Precautions.
2. Check that syringe tip cap is in place. Inspect the clarity of solution visually (Fig. 1a).
3. Depress plunger with tip cap to release the stopper seal (Fig. 1b). Unscrew tip cap from the syringe ensuring that there is no touch contamination of the syringe luer-connection (Fig. 2).
4. Push syringe plunger to expel the air. Do not use device if air is present.
5. After use, dispose of the syringe including any unused solution in accordance with institution's biohazardous waste policy or CDC guidelines (Fig. 6).

**WARNINGS AND PRECAUTIONS**  
0,9% Sodium Chloride solution is a sterile, aqueous solution having approximately the same osmotic pressure and composition as extracellular fluids. It is non-irritating to tissues and is used for vascular access devices in order to maintain catheter patency and to prevent contact between incompatible medications or fluids.

1. Read carefully instructions for use before using the product.
2. The product should be used only by qualified healthcare professionals.
3. Not to be used on Sterile Field.
4. Not to be used for injections as this is intended for flushing only.
5. Do not use product before use. Do not use if product is expired.
6. Do not use if tip cap/plunger is damaged, as it indicates syringe leakage or contamination.
7. Visually inspect the contents of each prefilled syringe for presence of saline solution, clear content, particulate matter, precipitate, discoloration and leakage prior to use. Do not use if any of the aforementioned conditions are present.
8. Do not allow air to be trapped in fluid path.
9. Do not re-sterilize. Discard the set after single use.
10. For single use only. The use of a single use syringe in multiple devices creates a potential risk for patient or user contamination and/or limited functionality of the device may lead to injury, illness or inflammation.

**STORAGE AND STABILITY**  
Prefilled Syringe with 0,9% sodium chloride solution should be stored between 5°C and 25°C, protect from direct sunlight and freezing.

**CAUTIONS**  
See Instructions for use

**REF** Product reference/ Art. No.

**LOT** Batch Number

**REF** For single use only

**LOT** Date of Manufacturing

**STERILE R** Use by/ Expiry date

**EU REPR** European Authorized Representative

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