Polyflush DS
0,9% sodium chlorid (0,9% NaCl)
Double sterile (DS) + Sterile Fluid + Sterile Syringe

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INSTRUCTIONS FOR USE

MATERIALS USED
Polypropylene (SBR), Sodium Chloride

INDICATION AND CLINICAL USE
Prefilled syringe with Saline solution is intended to be used only for syringing of medical vascular access devices (VADs) for maintaining patency of VADs.

CONTRAINDICATION
None Known

ADVERSE REACTION
Adverse reactions are known when the product is indicated.

DIRECTIONS FOR USE
Use aseptic technique throughout the procedure.

1. Open the pack and remove the syringe.
2. Visually inspect product for particulate matter and discoloration prior to use according to Warning and Precautions.
3. Check that syringe tip cap is in place. Inspect the clarity of solution visually (Fig. 1a).
4. Depress plunger with lip cap to replace the stopper seal (Fig. 1b).
5. Insert syringe tip cap in place. Inspect the clarity of solution visually (Fig. 1a).
6. Push syringe plunger to expel the air. Do not use device if no solution is observed when expelling the air (Fig. 3).
7. Connect Polyflush Syringe to vascular access device, taking care that there is no touch contamination of the connection (Fig. 2).
8. Push syringe plunger to fill the required volume of saline following instructions of policy (Fig. 5). In case of extreme plunger resistance, it is recommended that excess force is not exerted.
9. After use, dispose of the syringe including any unused solution in accordance with institution's biohazardous waste policy or CDC guidelines (Fig. 6).

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1. Åbn pakningen og fjern sprøjten.
2. Inspect visually the product for particulate matter and discoloration prior to use according to Warning and Precautions.
3. Check that syringe tip cap is in place. Inspect the clarity of solution visually (Fig. 1a).
4. Depress plunger with lip cap to replace the stopper seal (Fig. 1b).
5. Insert syringe tip cap in place. Inspect the clarity of solution visually (Fig. 1a).
6. Push syringe plunger to expel the air. Do not use device if no solution is observed when expelling the air (Fig. 3).
7. Connect Polyflush Syringe to vascular access device, taking care that there is no touch contamination of the connection (Fig. 2).
8. Push syringe plunger to fill the required volume of saline following instructions of policy (Fig. 5). In case of extreme plunger resistance, it is recommended that excess force is not exerted.
9. After use, dispose of the syringe including any unused solution in accordance with institution's biohazardous waste policy or CDC guidelines (Fig. 6).

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ACCIÓN
El 0,9% Solución Cloruro de Sodio es una solución salina having approximately the same osmotic pressure and composition as extracellular fluid. It is non-irritating to tissues. It is used to flush vascular access devices in order to maintain catheter patency and to prevent contact between incompatible medications or fluids.

WARNINGS AND PRECAUTIONS

1. Read the instructions for use before using the product.
2. The product should be used only by qualified healthcare professionals.
3. Do not use for injections as this is intended for flushing only.
4. Check the expiration date of product before use. Do not use if product is expired.
5. Do not use if tip cap/plunger is damaged, as it indicates syringe leakage or sterility.
6. Do not use if packaging or product has been damaged or contaminated.
7. Visually inspect the contents of each prefilled syringe for presence of salient, particulate matter, precipitate, discoloration and leakage prior to use. Do not use if any of the aforementioned is observed.
8. Do not use if seal path is damaged.
9. Do not re-sterilize. Discard the set after single use.
10. For single use in a single patient only. Re-use of single use devices creates a potential risk for patient or user. Contamination and/or limited functionality of the device may lead to injury, illness or infection for any direct incidental or consequential damage resulting from reuse of product.

STORAGE AND STABILITY
Polyflush Syringe with 0.9% sodium chloride solution should be stored between 5°C and 25°C, protect from direct sunlight and freezing.

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